

**BREATHE EASY,
SEE CLEARLY**

Build-in anti-fog strip
for clear vision

SafeMask
Antifog
Earloop Masks

AS 4381
LEVEL 3



A.R. Medicom (Australia) Pty Limited

5 Yulong Close, Moorebank NSW 2170

Tel: 1300 818 168

Email: cs@medicom-australia.com

Web: au.medicom.com

SafeMask

Antifog

Earloop Masks

Product Code

GMK222815

AS 4381
LEVEL 3

Packing

50 masks/box
10 boxes/case



Key Features and Benefits

FEATURES	BENEFITS
Adjustable, extra-length, enclosed aluminum nose piece with unique fog-free strip	Easily adjustable to form a strong seal across the nose for personalized fit; prevents fogging
High-quality filter layer meets AS4381:2015 Level 3 performance requirements	Optimal filtration for every situation or procedure
Non-woven, spun-bonded outer shell	Maximum fluid resistance without compromising breathability; will not lint, tear or shred
Ultrasonically sealed mask body	Eliminates rough edges that can irritate and scratch; ensures durable construction and eliminates holes and defects

Application



Medical



Dental



Aged
Care



Life
Science

Compliance

AS 4381:2015 Level 3 Standard	RESULTS
Bacterial filtration efficiency (BFE)	≥ 98%
Splash Resistant	@160 mm Hg
Differential Pressure (ΔP)	< 5.0 mm H ₂ O/cm ²

For more information contact:

Medicom Australia
 Tel: 1300 818 168 | E: cs@medicom-australia.com
au.medicom.com