


STERI*fit*[®]

Sterilisation Wrap



- SMS for strength and presents a tortuous pathway
 - Effective sterile barrier
- High resistance to tearing
 - High drapeability
- Available in 3 grades

STERIFIT[®] Sterilisation wrap is intended for wrapping medical devices that are to be sterilised and also to maintain sterility of the device until the packaging is opened for use, or until a predetermined shelf life is expired.

Dual Colour Standard Sterilisation Wrap

REORDER	DESCRIPTION	UNITS/CTN
STS001 [‡]	STERIFIT® Standard 30cm x 30cm	1000
STS002	STERIFIT® Standard 43cm x 43cm	500
STS003	STERIFIT® Standard 61cm x 61cm	500
STS004	STERIFIT® Standard 71cm x 71cm	300
STS005	STERIFIT® Standard 91.5cm x 91.5cm	150
STS006	STERIFIT® Standard 100cm x 106cm	150
STS007	STERIFIT® Standard 100cm x 140cm	100
STS008	STERIFIT® Standard 121.5cm x 121.5cm	100
STS009 [‡]	STERIFIT® Standard 130cm x 152.5cm	-
STS010	STERIFIT® Standard 137cm x 183cm	100

Dual Colour Heavy Duty Sterilisation Wrap

REORDER	DESCRIPTION	UNITS/CTN
STH001	Heavy Duty 71cm x 71cm	200
STH002	Heavy Duty 91.5cm x 91.5cm	150
STH003	Heavy Duty 106cm x 120cm	50
STH004 [‡]	Heavy Duty 106cm x 140cm	-
STH005	Heavy Duty 121.5cm x 121.5cm	50
STH006	Heavy Duty 130cm x 152.5cm	50
STH007	Heavy Duty 137cm x 183cm	50

Dual Colour Ultra Heavy Duty Sterilisation Wrap

REORDER	DESCRIPTION	UNITS/CTN
STU001	STERIFIT® Ultra Heavy Duty 106cm x 120cm	50
STU002 [‡]	STERIFIT® Ultra Heavy Duty 106cm x 140CM	-
STU004	STERIFIT® Ultra Heavy Duty 121.5cm x 121.5cm	50
STU005	STERIFIT® Ultra Heavy Duty 137cm x 183cm	40

[‡] Indent Only

Sterilisation Transportation Tray*

Sentry Medical Sterilisation Transportation Trays make transportation and storage of large or multiple sterile packs easier. The trays are used to prevent damage to the STERIFIT® wrapped packs before, during and after the sterilisation process. The trays can also be conveniently placed inside the autoclave during sterilisation, and are durable and reusable.

REORDER	DESCRIPTION	UNITS/CTN
STT001	560mm x 300mm x 35mm	20

* Available for loan through STERIFIT® agreements.

Contact customer service on 1300 995 999 or your local sales representative

Sentry Medical Pty Ltd
 22 Peter Brock Drive, Eastern Creek, NSW 2766, Australia
 Phone 1300 995 999 | Int Phone +612 8884 3995
 Fax 1300 995 998 | Int Fax +612 8884 3998
www.sentrymedical.com.au

